Hear now. And always

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The screenshot shows a web browser window displaying the website of the Norwegian Medical Association (Legeforeningen). The page is titled "Norsk forening for otorhinolaryngologi" and "DEN NORSKE LEGEFORENING". The main navigation menu includes "ORGANISASJON", "NETTUNDERVISNINGEN", "VEILEDER FOR ØRE-NESE-HALSFAGET", and "MØTER". The current page is "Audiologi" under the "Veileder for øre-nese-halsfaget / Audiologi" section. The content includes the title "Audiologi", a subtitle "Audiologisk utvalg 2011, ved Haakon Arnesen, Niels Chr. Stenklev, Geir Siem. Utvalg for kvalitetssikring 2011, ved redaktør Veileder for Øre-nese-halsfaget Torgrim Fuhr", and the date "08. mars 2012. Emner: øre-nese-halssykdommer". There are also social media sharing options for "Anbefal" and "Bli den første av dine venner til å anbefale dette". On the right side, there is a "Les også" section with a list of related articles.

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
Hear now. And always

# Hvor godt fungerer personer med CI?

# Assessment of Speech Understanding After Cochlear Implantation in Adult Hearing Aid Users

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Buchman C et al. Assessment of Speech Understanding after Cochlear Implantation in Adult Hearing Aid Users. JAMA Otolaryngol- Head Neck Surg. Published online Aug 27th 2020

© Cochlear Limited 2023 D2072935-V1

Research

JAMA Otolaryngology-Head & Neck Surgery | Original Investigation

## Assessment of Speech Understanding After Cochlear Implantation in Adult Hearing Aid Users A Nonrandomized Controlled Trial

Craig A. Buchman, MD; Jacques A. Herzog, MD; Jonathan L. McLunkin, MD; Cameron C. Wick, MD; Nedim Durakovic, MD; Jill B. Firsz, PhD; Dorina Kallogjeri, MD, MPH; for the CIS32 Study Group

**IMPORTANCE** Cochlear implants were approved for use in adults in the 1980s, but use remains low owing to a lack of awareness regarding cochlear implantation candidacy criteria and expected outcomes. There have been limited, small series examining the safety and effectiveness of cochlear implantation in adult hearing aid (HA) users with and without mild cognitive impairment (MCI).

**OBJECTIVE** To investigate the safety and effectiveness of a single-ear cochlear implant in a group of optimized adult HA users with and without MCI across a variety of domains.

**DESIGN, SETTING, AND PARTICIPANTS** In this nonrandomized controlled trial, a multicenter, prospective, repeated-measures investigation was conducted at 13 US institutions. The setting was academic and community-based cochlear implant programs. Eligible participants were 100 adults (aged >18 years) with postlinguistic onset of bilateral moderate sloping to profound or worse sensorineural hearing loss ( $\leq 20$  years' duration). Fluent English speakers underwent an optimized bilateral HA trial for at least 30 days. Individuals with aided Consonant-Vowel Nucleus-Consonant (CNC) word score in quiet of 40% or less correct in the ear to be implanted and 50% or less correct in the contralateral ear were offered cochlear implants. The first participant was enrolled on February 20, 2017, and the last participant was enrolled on May 3, 2018. The final follow-up was on December 21, 2018.

**INTERVENTIONS** Participants received the same cochlear implant system and contralateral HA.

**MAIN OUTCOMES AND MEASURES** The primary outcome measure was speech understanding in quiet (CNC word score) using both the cochlear implant and opposite ear HA. Secondary outcome measures included the following: adverse events; speech understanding in noise (AzBio signal-to-noise ratio of +10 dB [+10 SNR]) Health Utilities Index Mark 3 (HUI3); Speech, Spatial, and Qualities of Hearing Questionnaire 49 (SSQ49); and Montreal Cognitive Assessment (MoCA).

**RESULTS** The median age at cochlear implantation of the 96 patients included in the trial was 71 years (range, 23-91 years), and 62 patients (65%) were male. Three serious adverse events requiring revision surgery occurred, and all resolved without sequelae. By 6 months after activation, the absolute marginal mean change in CNC word score and AzBio +10 SNR was 40.5% (95% CI, 35.9%-45.0%) and 24.1% (95% CI, 18.9%-29.4%), respectively. Ninety-one percent (87 of 96) of participants had a clinically important improvement ( $>15\%$ ) in the CNC word score in the implant ear. Mild cognitive impairment (MoCA total score  $\leq 25$ ) was observed in 48 of 81 study participants (59%) at baseline. Speech perception marginal mean improvements were similar between individuals with and without baseline MCI, with values of 40.9% (95% CI, 35.2%-46.6%) and 39.6% (95% CI, 31.8%-47.4%), respectively, for CNC word score and 27.5% (95% CI, 21.0%-33.9%) and 17.8% (95% CI, 9.0%-26.6%), respectively, for AzBio +10 SNR. Statistically significant and clinically important improvements in the HUI3 and SSQ49 were evident at 6 months.

**CONCLUSIONS AND RELEVANCE** The findings of this nonrandomized controlled trial seem to indicate that cochlear implants are safe and effective in restoring speech understanding in both quiet and noise and improve quality of life in individuals with and without MCI.

**TRIAL REGISTRATION** ClinicalTrials.gov Identifier: NCT03007472

JAMA Otolaryngol Head Neck Surg. doi:10.1001/jamaoto.2020.1584  
Published online August 27, 2020.

- Viewpoint and Invited Commentary
- Author Audio Interview
- Related articles
- Supplemental content

**Author Affiliations:** Department of Otolaryngology-Head & Neck Surgery, Washington University in St. Louis, St. Louis, Missouri (Buchman, Herzog, McLunkin, Wick, Durakovic, Firsz, Kallogjeri); Statistics Editor, JAMA Otolaryngology-Head & Neck Surgery (Kallogjeri).

**Group Information:** The members of the CIS32 Study Group appear at the end of the article.

**Corresponding Author:** Craig A. Buchman, MD, Department of Otolaryngology-Head & Neck Surgery, Washington University School of Medicine in St. Louis, 660 S Euclid Ave, Campus Box 8115, St. Louis, MO 63110 (buchman@wustl.edu).

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
# Hearing and Quality-of-Life Outcomes After Cochlear Implantation in Adult Hearing Aid Users 65 Years or Older

A Secondary Analysis of a Nonrandomized Clinical Trial

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Wick C et al. Hearing and QoL outcomes after Cochlear Implantation in Adult Hearing Aid users 65 years or older. JAMA Otolaryngol Head Neck Surg Published online Aug 27th 2020

## Hearing and Quality-of-Life Outcomes After Cochlear Implantation in Adult Hearing Aid Users 65 Years or Older A Secondary Analysis of a Nonrandomized Clinical Trial

Cameron C. Wick, MD; Dorina Kallgjeri, MD; Jonathan L. McJunkin, MD; Nedim Durakovic, MD; Laura K. Holden, AuD; Jacques A. Herzog, MD; Jill B. Firszt, PhD; Craig A. Buchman, MD; for the CIS32 Study Group

**IMPORTANCE** Hearing loss, especially moderate to severe forms, has the potential to negatively affect an individual's physical, social, emotional, and cognitive well-being. Moreover, having ineffective binaural hearing increases difficulty understanding speech in noise and leads to a greater degree of social isolation and loneliness and a reduced quality of life (QoL).

**OBJECTIVE** To explore the audiometric and holistic effects of cochlear implantation in a group of adults 65 years or older compared with an optimized bilateral hearing aid condition.

**DESIGN, SETTING, AND PARTICIPANTS** This ad hoc secondary analysis of a prospective, single-subject, repeated-measures nonrandomized clinical trial included 13 cochlear implantation centers across the United States. Participants 65 years or older with postlingual bilateral moderate-to-profound sensorineural hearing loss with aided Consonant-Vowel Nucleus-Consonant (CNC) word scores in quiet of 40% or less in the ear to undergo implantation and 50% or less in the contralateral ear were included in the analysis. Baseline QoL testing was performed after 1 month of optimized bilateral hearing aid use. Participants were enrolled from February 20, 2017, to May 3, 2018, and follow-up was completed December 21, 2018. Data were analyzed from March 25, 2019, to March 31, 2020.

**INTERVENTIONS** Unilateral implantation with a slim, modular cochlear implant device. Hearing aid use in the contralateral ear was required through the 6-month primary end-point interval.

**MAIN OUTCOMES AND MEASURES** The primary objective was to evaluate speech perception before and 6 months after activation of a new cochlear implant. Secondary objectives were QoL metrics in the everyday listening condition before and 6 months after implantation.

**RESULTS** Seventy participants (51 men [73%]) with a median age of 74 (range, 65–91) years were included in the analysis. No major adverse events occurred. Mixed-model analysis with estimated marginal means and 95% CIs compared preimplantation baseline performance with 6-month postimplantation performance. A clinically important improvement in CNC words was shown in the bimodal condition, with a mean difference of 37.2% (95% CI, 32.0%–42.4%), and in the unilateral (cochlear implant only) condition, with a mean difference of 44.3% (95% CI, 39.0%–49.2%). A clinically important improvement in noise (AzBio sentences signal-to-noise ratio of +10 dB) was also shown, with a mean difference of 21.6% (95% CI, 15.7%–27.5%) in the bimodal condition and 24.5% (95% CI, 18.3%–30.7%) in the unilateral condition. The Health Utilities Index Mark 3 multiple-attribute score improved by 0.186 (95% CI, 0.136–0.234); the Speech, Spatial, and Qualities of Hearing Scale total score improved by 2.58 (95% CI, 2.18–2.99); and a novel Device Use Questionnaire reported 94% of participants were satisfied with overall hearing in the everyday listening condition.

**CONCLUSIONS AND RELEVANCE** This subgroup analysis of patients 65 years or older enrolled in a within-subject clinical trial of cochlear implantation demonstrated clinically meaningful audiometric and QoL benefit with an acceptable risk profile. These findings suggest that cochlear implantation in older adults may facilitate the concept of healthy aging.

**TRIAL REGISTRATION** ClinicalTrials.gov Identifier: NCT03007472

Viewpoint and Invited Commentary

Author Audio Interview

Related articles

Supplemental content

**Author Affiliations:** Department of Otolaryngology–Head and Neck Surgery, Washington University School of Medicine in St. Louis, St. Louis, Missouri (Wick, Kallgjeri, McJunkin, Durakovic, Holden, Herzog, Firszt, Buchman); Statistics Editor, JAMA Otolaryngology–Head & Neck Surgery (Kallgjeri).

**Group Information:** A complete list of the members of the CIS32 Study Group and participating centers appears at the end of this article.


**Corresponding Author:** Cameron C. Wick, MD, Department of Otolaryngology–Head and Neck Surgery, Washington University School of Medicine in St. Louis, 660 S Euclid Ave, PO Box 8115, St. Louis, MO 63110 (cameron.wick@wustl.edu).

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
2. Buchman et al. International consensus paper: Unilateral Cochlear Implants for Severe, Profound, or Moderate Sloping to Profound Sensorineural Hearing Loss: A Systematic Review and Consensus Statements. *JAMA Otolaryngol Head Neck Surg.*

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... kan være sosial

... får bedre livskvalitet

\* 1.4 The impact of unaddressed hearing loss, WHO's World report on hearing, ISBN 978-92-4-002048-1

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